

# SEJONG TABLET PRESS REPORT

(1/2)

Software version ( 2.0 )

## PRODUCTION REPORT

### Product Information

|                  |                  |                      |                      |
|------------------|------------------|----------------------|----------------------|
| Station No :     | Station 1        |                      |                      |
| Machine Name :   | MC081 SEJONG 49D | Product Name :       | Mirtazapine OD 15mg  |
| User ID :        | pk8943           | Batch NO. :          | MPA26053             |
| Print Interval : | 385000 Tabs      | Running Time :       | 02 Hour 03 Min45 Sec |
| Total Counter :  | Tabs             | Total Running Time : | 270 Hour 40 Min      |

### Setting Value

|                          |                   |                        |                       |
|--------------------------|-------------------|------------------------|-----------------------|
| Feeder :                 | Auto ; 40 %       | Manual ; 9 RPM         |                       |
| Filling Cam :            | Cam B 4.5-10mm    | Target Quantity :      | 1540000 Tabs          |
| Air Pressure Low Limit : | 400 Kpa           | Hydraulic High Limit : | 15 Mpa (Low: 0.1 Mpa) |
| Oil Lubrication S1 :     | Interval: 300 Min | Supply: 3 Sec          | Upper Punch           |
| Oil Lubrication S2 :     | Interval: 300 Min | Supply: 3 Sec          | Lower Punch           |
| Oil Lubrication S3 :     | Interval: 999 Min | Supply: 0 Sec          | Lower Head            |
| Powder Supply Time :     | 30 Sec            | Initial Reject Time :  | 6 Sec                 |

### Control Limits & Stop Conditions

| Parameter                   | % Setting | kN Limit | Stop Condition   |
|-----------------------------|-----------|----------|------------------|
| HSP (High Stop Pressure)    | 45 %      | 8.8 kN   | Stop: Yes        |
| HEP (High Error Pressure)   | 28 %      | 7.77 kN  | 10 Rot / 10 Tabs |
| HCP (High Control Pressure) | 3 %       | 6.25 kN  | 5 Times          |
| Ref (Reference Pressure)    | -         | 6.07 kN  | -                |
| LCP (Low Control Pressure)  | 3 %       | 5.89 kN  | 5 Times          |
| LEP (Low Error Pressure)    | 27 %      | 4.43 kN  | 10 Rot / 10 Tabs |
| LSP (Low Stop Pressure)     | 45 %      | 3.34 kN  | Stop: Yes        |
| SD Limit                    | 15 %      | -        | Stop: Yes        |
| Pre HSP                     | -         | 9.91 kN  | Stop: Yes        |

### Pressure Data

| Section                   | Mean (kN)          | SD (%) | Min (kN) [Punch] | Max (kN) [Punch] |
|---------------------------|--------------------|--------|------------------|------------------|
| Pre Pressure              | 2.38 kN            | 0.7 %  | 2.32 kN [#32]    | 2.44 kN [#13]    |
| Main Pressure             | 6 kN               | 2.1 %  | 5.63 kN [#62]    | 6.49 kN [#3]     |
| Filling Depth Adjustments | Increase: 11 times |        |                  |                  |

Date: 2026-09-13 00:00:00

Operator: p.s.krishna murthy

Signature: \_\_\_\_\_

# SEJONG TABLET PRESS REPORT

(2/2)

Software version ( 2.0 )

## PRODUCTION REPORT - OPERATIONS & TABLET DATA

### Operation Value

|                           |          |                        |                          |
|---------------------------|----------|------------------------|--------------------------|
| Disk Speed :              | 21 RPM   | Capacity :             | 183960 Tabs/hour         |
| Feeder Status :           | AUTO     | Feeder Speed :         | 8 RPM                    |
| Pre-Pressure Thickness :  | 11.23 mm | Lower Punch Position : | 13.23 mm (Penet: 2 mm)   |
| Main-Pressure Thickness : | 4.19 mm  | Lower Punch Position : | 8.19 mm (Penet: 4 mm)    |
| Filling Depth :           | 7.69 mm  | Current Cam :          | Cam B 4.5-10mm           |
| Main Air Pressure :       | 661 Kpa  | Hydraulic Pressure :   | 7.6 Mpa                  |
| Oil S1 Remain Time :      | 56 Min   | Oil S2 / S3 Remain :   | S2: 56 Min   S3: 653 Min |
| Powder Status :           | Enable   | Dust Collector :       | ON                       |
| Initial Reject :          | ON       | Buzzer :               | ON                       |

### Tablet Data & Production Counters

| Metric                        | Value / Count (Tabs) | Percentage / Details        |
|-------------------------------|----------------------|-----------------------------|
| Total Counter                 | 385,158 Tabs         | 100% Machine Processed      |
| A.W.C. Counter                | 366,022 Tabs         | Automatic Weight Controlled |
| HEP (High Error Punch Reject) | 0 Tabs               | 0 (0%)                      |
| LEP (Low Error Punch Reject)  | 0 Tabs               | 0 (0%)                      |
| GOOD Tablets                  | 366,022 Tabs         | 366022 (100.0%)             |

Date: 2026-09-13 00:00:00

Operator: p.s.krishna murthy

Signature: \_\_\_\_\_