

SEJONG TABLET PRESS REPORT

(1/2)

Software version (2.0)

PRODUCTION REPORT

Product Information

Station No :	Station 1		
Machine Name :	MC081 SEJONG 49D	Product Name :	Mirtazapine OD 15mg
User ID :	pk8943	Batch NO. :	MPA26060
Print Interval :	385000 Tabs	Running Time :	01 Hour 58 Min34 Sec
Total Counter :	Tabs	Total Running Time :	326 Hour 23 Min

Setting Value

Feeder :	Auto ; 40 %	Manual ; 9 RPM	
Filling Cam :	Cam B 4.5-10mm	Target Quantity :	1540000 Tabs
Air Pressure Low Limit :	400 Kpa	Hydraulic High Limit :	15 Mpa (Low: 0.1 Mpa)
Oil Lubrication S1 :	Interval: 300 Min	Supply: 3 Sec	Upper Punch
Oil Lubrication S2 :	Interval: 300 Min	Supply: 3 Sec	Lower Punch
Oil Lubrication S3 :	Interval: 999 Min	Supply: 0 Sec	Lower Head
Powder Supply Time :	30 Sec	Initial Reject Time :	5.8 Sec

Control Limits & Stop Conditions

Parameter	% Setting	kN Limit	Stop Condition
HSP (High Stop Pressure)	45 %	8.36 kN	Stop: Yes
HEP (High Error Pressure)	28 %	7.38 kN	10 Rot / 10 Tabs
HCP (High Control Pressure)	3 %	5.94 kN	5 Times
Ref (Reference Pressure)	-	5.77 kN	-
LCP (Low Control Pressure)	3 %	5.6 kN	5 Times
LEP (Low Error Pressure)	27 %	4.21 kN	10 Rot / 10 Tabs
LSP (Low Stop Pressure)	45 %	3.17 kN	Stop: Yes
SD Limit	15 %	-	Stop: Yes
Pre HSP	-	9.95 kN	Stop: Yes

Pressure Data

Section	Mean (kN)	SD (%)	Min (kN) [Punch]	Max (kN) [Punch]
Pre Pressure	3.81 kN	0.6 %	3.75 kN [#72]	3.88 kN [#28]
Main Pressure	5.84 kN	1.8 %	5.59 kN [#36]	6.11 kN [#21]
Filling Depth Adjustments	Increase: 0 times			

Date: 2026-09-20 00:00:00

Operator: p.s.krishna murthy

Signature: _____

SEJONG TABLET PRESS REPORT

(2/2)

Software version (2.0)

PRODUCTION REPORT - OPERATIONS & TABLET DATA

Operation Value

Disk Speed :	22 RPM	Capacity :	192720 Tabs/hour
Feeder Status :	AUTO	Feeder Speed :	8 RPM
Pre-Pressure Thickness :	7.22 mm	Lower Punch Position :	9.22 mm (Penet: 2 mm)
Main-Pressure Thickness :	4.18 mm	Lower Punch Position :	8.18 mm (Penet: 4 mm)
Filling Depth :	7.92 mm	Current Cam :	Cam B 4.5-10mm
Main Air Pressure :	697 Kpa	Hydraulic Pressure :	7.6 Mpa
Oil S1 Remain Time :	21 Min	Oil S2 / S3 Remain :	S2: 21 Min S3: 315 Min
Powder Status :	Enable	Dust Collector :	ON
Initial Reject :	ON	Buzzer :	ON

Tablet Data & Production Counters

Metric	Value / Count (Tabs)	Percentage / Details
Total Counter	385,052 Tabs	100% Machine Processed
A.W.C. Counter	373,833 Tabs	Automatic Weight Controlled
HEP (High Error Punch Reject)	0 Tabs	0 (0%)
LEP (Low Error Punch Reject)	0 Tabs	0 (0%)
GOOD Tablets	373,833 Tabs	373833 (100.0%)

Date: 2026-09-20 00:00:00

Operator: p.s.krishna murthy

Signature: _____