

# SEJONG TABLET PRESS REPORT

(1/2)

Software version ( 2.0 )

## PRODUCTION REPORT

### Product Information

|                  |                  |                      |                      |
|------------------|------------------|----------------------|----------------------|
| Station No :     | Station 1        |                      |                      |
| Machine Name :   | MC081 SEJONG 49D | Product Name :       | Mirtazapine OD 15mg  |
| User ID :        | pk8943           | Batch NO. :          | MPA26054             |
| Print Interval : | 385000 Tabs      | Running Time :       | 06 Hour 11 Min12 Sec |
| Total Counter :  | Tabs             | Total Running Time : | 282 Hour 57 Min      |

### Setting Value

|                          |                   |                        |                       |
|--------------------------|-------------------|------------------------|-----------------------|
| Feeder :                 | Auto ; 40 %       | Manual ; 9 RPM         |                       |
| Filling Cam :            | Cam B 4.5-10mm    | Target Quantity :      | 1540000 Tabs          |
| Air Pressure Low Limit : | 400 Kpa           | Hydraulic High Limit : | 15 Mpa (Low: 0.1 Mpa) |
| Oil Lubrication S1 :     | Interval: 300 Min | Supply: 3 Sec          | Upper Punch           |
| Oil Lubrication S2 :     | Interval: 300 Min | Supply: 3 Sec          | Lower Punch           |
| Oil Lubrication S3 :     | Interval: 999 Min | Supply: 0 Sec          | Lower Head            |
| Powder Supply Time :     | 30 Sec            | Initial Reject Time :  | 6 Sec                 |

### Control Limits & Stop Conditions

| Parameter                   | % Setting | kN Limit | Stop Condition   |
|-----------------------------|-----------|----------|------------------|
| HSP (High Stop Pressure)    | 45 %      | 9.56 kN  | Stop: Yes        |
| HEP (High Error Pressure)   | 28 %      | 8.44 kN  | 10 Rot / 10 Tabs |
| HCP (High Control Pressure) | 3 %       | 6.79 kN  | 5 Times          |
| Ref (Reference Pressure)    | -         | 6.59 kN  | -                |
| LCP (Low Control Pressure)  | 3 %       | 6.39 kN  | 5 Times          |
| LEP (Low Error Pressure)    | 27 %      | 4.81 kN  | 10 Rot / 10 Tabs |
| LSP (Low Stop Pressure)     | 45 %      | 3.63 kN  | Stop: Yes        |
| SD Limit                    | 15 %      | -        | Stop: Yes        |
| Pre HSP                     | -         | 9.91 kN  | Stop: Yes        |

### Pressure Data

| Section                   | Mean (kN)          | SD (%) | Min (kN) [Punch] | Max (kN) [Punch] |
|---------------------------|--------------------|--------|------------------|------------------|
| Pre Pressure              | 2.45 kN            | 0.8 %  | 2.39 kN [#61]    | 2.51 kN [#5]     |
| Main Pressure             | 6.7 kN             | 2.9 %  | 6.28 kN [#29]    | 7.4 kN [#18]     |
| Filling Depth Adjustments | Increase: 15 times |        |                  |                  |

Date: 2026-09-15 00:00:00

Operator: p.s.krishna murthy

Signature: \_\_\_\_\_

# SEJONG TABLET PRESS REPORT

(2/2)

Software version ( 2.0 )

## PRODUCTION REPORT - OPERATIONS & TABLET DATA

### Operation Value

|                           |          |                        |                           |
|---------------------------|----------|------------------------|---------------------------|
| Disk Speed :              | 21 RPM   | Capacity :             | 183960 Tabs/hour          |
| Feeder Status :           | AUTO     | Feeder Speed :         | 8 RPM                     |
| Pre-Pressure Thickness :  | 11.23 mm | Lower Punch Position : | 13.23 mm (Penet: 2 mm)    |
| Main-Pressure Thickness : | 4.12 mm  | Lower Punch Position : | 8.12 mm (Penet: 4 mm)     |
| Filling Depth :           | 7.52 mm  | Current Cam :          | Cam B 4.5-10mm            |
| Main Air Pressure :       | 579 Kpa  | Hydraulic Pressure :   | 7.5 Mpa                   |
| Oil S1 Remain Time :      | 219 Min  | Oil S2 / S3 Remain :   | S2: 219 Min   S3: 915 Min |
| Powder Status :           | Enable   | Dust Collector :       | ON                        |
| Initial Reject :          | ON       | Buzzer :               | ON                        |

### Tablet Data & Production Counters

| Metric                        | Value / Count (Tabs) | Percentage / Details        |
|-------------------------------|----------------------|-----------------------------|
| Total Counter                 | 1,155,188 Tabs       | 100% Machine Processed      |
| A.W.C. Counter                | 1,096,168 Tabs       | Automatic Weight Controlled |
| HEP (High Error Punch Reject) | 0 Tabs               | 0 (0%)                      |
| LEP (Low Error Punch Reject)  | 0 Tabs               | 0 (0%)                      |
| GOOD Tablets                  | 1,096,168 Tabs       | 1096168 (100.0%)            |

Date: 2026-09-15 00:00:00

Operator: p.s.krishna murthy

Signature: \_\_\_\_\_